

Product Datasheet

Phencyclidine (PCP)-BSA Antigen LA502

Product Name	Phencyclidine (PCP)-BSA Antigen
Catalog Number	LA502
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	$\geq 90\%$ by SDS-PAGE
Buffer	10 mM Phosphate Buffer, pH 7.2
Matched Pairs	HM927 (detection)
Preservatives	0.01% Sodium Azide