

Product Datasheet

Carcinoembryonic Antigen (CEA), Native >95% LA505L

Product Name	Carcinoembryonic Antigen (CEA), Native >95%
Catalog Number	LA505L
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Lyophilized Native Antigen
Storage	Lyophilized or Reconstituted: Short Term (<1 week): 2-8°C Long Term: -20°C. Avoid repeated freezing and thawing
Purity	>95% by SDS-PAGE
Source	Human Liver
Molecular Weight	180 kDa