

Product Datasheet

Respiratory Syncytial Virus (RSV-fusion) Antigen, Recombinant >90% LA512

| | |
|----------------------------|---|
| Product Name | Respiratory Syncytial Virus (RSV-fusion) Antigen, Recombinant >90% |
| Catalog Number | LA512 |
| Tested Applications | ELISA, LF |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Recombinant Antigen |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Purity | >90% |
| Buffer | Acetate Buffer with salts and other additives, pH 5.0 |
| Matched Pairs | HM947 |
| Expression System | Mammalian Cell |
| Preservatives | 0.09% Sodium Azide |