

## Product Datasheet

### Respiratory Syncytial Virus (RSV) Antigen, Native >90% LA513

<b>Product Name</b>	Respiratory Syncytial Virus (RSV) Antigen, Native >90%
<b>Catalog Number</b>	LA513
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Inactivated Native Antigen
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Source</b>	Hep-2 cells
<b>Buffer</b>	10 mM Phosphate Buffer with additives and salts, pH 7.4 ± 0.2
<b>Matched Pairs</b>	HM947
<b>Preservatives</b>	0.09% Sodium Azide