

Product Datasheet

Respiratory Syncytial Virus (RSV) Antigen, Native >90% LA513

Product Name	Respiratory Syncytial Virus (RSV) Antigen, Native >90%
Catalog Number	LA513
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Inactivated Native Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	Hep-2 cells
Buffer	10 mM Phosphate Buffer with additives and salts, pH 7.4 ± 0.2
Matched Pairs	HM947
Preservatives	0.09% Sodium Azide