

Product Datasheet

Adenovirus Hexon Antigen, Recombinant $\geq 90\%$ LA514

Product Name	Adenovirus Hexon Antigen, Recombinant $\geq 90\%$
Catalog Number	LA514
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	$\geq 90\%$
Buffer	Tris buffer with additives and salts, pH 8.0 \pm 0.2
Matched Pairs	HM948
Expression System	Mammalian Cell
Preservatives	0.09% Sodium Azide