

Product Datasheet

Adenovirus Hexon Antigen, Recombinant $\geq 90\%$ LA514

| | |
|----------------------------|---|
| Product Name | Adenovirus Hexon Antigen, Recombinant $\geq 90\%$ |
| Catalog Number | LA514 |
| Tested Applications | ELISA, LF |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Recombinant Antigen |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Purity | $\geq 90\%$ |
| Buffer | Tris buffer with additives and salts, pH 8.0 \pm 0.2 |
| Matched Pairs | HM948 |
| Expression System | Mammalian Cell |
| Preservatives | 0.09% Sodium Azide |