

Product Datasheet

Escherichia coli (E. coli) O157 Verotoxin 1 Antigen, Recombinant >90% LA519

Product Name	Escherichia coli (E. coli) O157 Verotoxin 1 Antigen, Recombinant >90%
Catalog Number	LA519
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	Carbonate Buffer, pH 8.3
Expression System	E. coli
Preservatives	0.09% Sodium Azide