

## Product Datasheet

### Human Papilloma Virus (HPV-16) Antigen, Recombinant >90% LA523

<b>Product Name</b>	Human Papilloma Virus (HPV-16) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA523
<b>Tested Applications</b>	Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (GST-tag)
<b>Storage</b>	Short Term ( $\leq$ 1 week) 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>90% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, 100 mM Arginine and 3M Urea
<b>Molecular Weight</b>	83 kDa
<b>Expression System</b>	E. coli