

# Product Datasheet

## Human Perforin ELISA Set 851.860.005

Product Name	Human Perforin ELISA Set
Catalog Number	851.860.005
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	62.5-2000 pg/ml
Sensitivity	40 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human Perforin
Sample Types	Serum, Plasma, Cell culture supernatant