

# Product Datasheet

## Cholera Toxin Subunit B Antigen, Recombinant >90% LA526L

Product Name	Cholera Toxin Subunit B Antigen, Recombinant >90%
Catalog Number	LA526L
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Lyophilized Recombinant Antigen (His-tag)
Storage	Lyophilized: Short Term 2-8°C Long Term: -20°C. Reconstituted: -20°C. Avoid repeated freezing and thawing
Note	Reconstitute with sterile water with gentle mixing. Addition of 0.1% BSA may increase stability.
Purity	>90% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Expression System	Mammalian Cell