

# Product Datasheet

## Diphtheria Toxin Antigen, Native >98% LA527L

Product Name	Diphtheria Toxin Antigen, Native >98%
Catalog Number	LA527L
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Lyophilized Native Antigen
Storage	Lyophilized: 2-8°C. Reconstituted: Short Term (≤ 2 weeks) 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
Note	Reconstitute with distilled water.
Purity	>98%
Source	Corynebacterium diphtheria culture
Molecular Weight	58 kDa