

Product Datasheet

Human Perforin ELISA Set 851.860.010

Product Name	Human Perforin ELISA Set
Catalog Number	851.860.010
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	62.5-2000 pg/ml
Sensitivity	40 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human Perforin
Sample Types	Serum, Plasma, Cell culture supernatant