

Product Datasheet

Heparin Binding Protein (HBP) Antigen, Recombinant LA537

Product Name Heparin Binding Protein (HBP) Antigen, Recombinant

Catalog Number LA537

Tested Applications Control

Alternative Names Azurocidin

Brand EastCoast Bio

Form/Appearance Purified Human Recombinant Antigen

Storage -20°C

Purity >95% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Matched Pairs HM721, HM722

Expression System Mammalian Cell

Preservatives 0.09% Sodium Azide