

Product Datasheet

Barbiturate-BSA Antigen LA538

Product Name	Barbiturate-BSA Antigen
Catalog Number	LA538
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservatives	0.02% Sodium Azide