

Product Datasheet

Cotinine-BSA Antigen LA539

| | |
|----------------------------|---|
| Product Name | Cotinine-BSA Antigen |
| Catalog Number | LA539 |
| Tested Applications | ELISA, LF |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid; Antigen-BSA conjugate |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.2 |
| Matched Pairs | PP195IF (detection) |
| Preservatives | 0.02% Sodium Azide |