

Product Datasheet

Cotinine-BSA Antigen LA539

Product Name	Cotinine-BSA Antigen
Catalog Number	LA539
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	PP195IF (detection)
Preservatives	0.02% Sodium Azide