

Product Datasheet

Estrone-3-Sulfate (E1S)-BSA Antigen LA551

Product Name	Estrone-3-Sulfate (E1S)-BSA Antigen
Catalog Number	LA551
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	≥90% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM1016 (detection)
Preservatives	0.01% Sodium Azide