

# Product Datasheet

## Human IP-10 ELISA Set 851.870.010

Product Name	Human IP-10 ELISA Set
Catalog Number	851.870.010
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	6.25 -200 pg/ml
Sensitivity	5.7 pg/ml
Tested Applications	ELISA
Alternative Names	CXCL10
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human molecules
Specificity	Recognizes both natural and recombinant form of the molecule
Sample Types	Serum, Plasma, Cell culture supernatant, buffered solutions