

## Product Datasheet

### Hepatitis Be Antigen (HBeAg VLP), Recombinant >95% LA577

<b>Product Name</b>	Hepatitis Be Antigen (HBeAg VLP), Recombinant >95%
<b>Catalog Number</b>	LA577
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Recombinant Virus Like Particle Antigen, (His-tag)
<b>Storage</b>	Short Term ( $\leq 1$ week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Molecular Weight</b>	18 kDa
<b>Expression System</b>	E. coli
<b>Preservatives</b>	0.05% Sodium Azide