

Product Datasheet

Anti-Mullerian Hormone (AMH 19-449) Antigen, Recombinant >95% LA578

Product Name	Anti-Mullerian Hormone (AMH 19-449) Antigen, Recombinant >95%
Catalog Number	LA578
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 19-449
Purity	>95% by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.0 ± 0.2
Molecular Weight	50 kDa
Matched Pairs	HM1050
Expression System	Mammalian Cell
Preservatives	0.09% Sodium Azide