

## Product Datasheet

### Cannabinoid (THC)-BSA Antigen LA581

<b>Product Name</b>	Cannabinoid (THC)-BSA Antigen
<b>Catalog Number</b>	LA581
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥90% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 ± 0.2
<b>Matched Pairs</b>	HM1062 (detection)
<b>Preservatives</b>	0.01% Sodium Azide