

Product Datasheet

Aripiprazole-BSA Antigen LA799

Product Name	Aripiprazole-BSA Antigen
Catalog Number	LA799
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM1262 or HM1263 (detection)
Preservatives	0.05% Sodium Azide