

Product Datasheet

Coronavirus (COVID-19 Spike Protein S1; 1-681) Antigen, Recombinant >95% LA610

Product Name	Coronavirus (COVID-19 Spike Protein S1; 1-681) Antigen, Recombinant >95%
Catalog Number	LA610
Tested Applications	ELISA, LF, Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 1-681 of spike protein S1 domain.
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Molecular Weight	76 kDa
Expression System	Mammalian Cell (HEK293)
Preservatives	0.09% Sodium Azide