

Product Datasheet

Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90% LA615

Product Name	Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90%
Catalog Number	LA615
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term (<2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing
Note	Peptide composed of amino acids 1-419 of nucleoprotein.
Purity	>90%
Buffer	25 mM Tris-HCl, pH 8.0
Molecular Weight	48 kDa
Expression System	E. coli