

## Product Datasheet

### Coronavirus (COVID-19 Spike Protein; Full Length) Antigen, Recombinant >90% LA636

<b>Product Name</b>	Coronavirus (COVID-19 Spike Protein; Full Length) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA636
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Full Length Recombinant Antigen
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Peptide composed of full length spike protein sequence.
<b>Purity</b>	≥90%
<b>Buffer</b>	Tris Buffer with salts and additives, pH 8.0
<b>Molecular Weight</b>	138 kDa
<b>Expression System</b>	Mammalian Cell (HEK293)