

Product Datasheet

Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90% LA637

Product Name Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90%

Catalog Number LA637

Tested Applications ELISA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Note Peptide composed of amino acids 1-419 of nucleoprotein.

Purity $\geq 90\%$

Buffer Tris buffer with salts and additives, pH 8.0

Molecular Weight 46.7 kDa

Expression System Mammalian Cell (HEK293)

Preservatives 0.09% Sodium Azide