

Product Datasheet

Human IL-27 ELISA Set 851.890.001

| | |
|---------------------|--|
| Product Name | Human IL-27 ELISA Set |
| Catalog Number | 851.890.001 |
| Description | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range | 31.3-1000 pg/ml |
| Sensitivity | 12 pg/ml |
| Tested Applications | ELISA |
| Alternative Names | NULL |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human cytokines |
| Specificity | Recognizes both natural and recombinant human IL-27 |
| Sample Types | Supernatants, buffered solutions, serum, plasmas and other body fluids. |