

Product Datasheet

Carisoprodol (Soma)-BSA Antigen LA648

Product Name	Carisoprodol (Soma)-BSA Antigen
Catalog Number	LA648
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	P02-99-72R-IF or P02-99-73R-IF (detection)
Preservatives	0.02% Sodium Azide