

Product Datasheet

Angiotensin-converting Enzyme 2 (ACE2) Antigen, Recombinant $\geq 90\%$ LA650

Product Name	Angiotensin-converting Enzyme 2 (ACE2) Antigen, Recombinant $\geq 90\%$
Catalog Number	LA650
Tested Applications	ELISA, LF, CMIA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (hFc-tag)
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	$\geq 90\%$
Buffer	10 mM Phosphate Buffer, 150 mM Sodium Chloride, pH 7.4
Molecular Weight	250 kDa
Expression System	Mammalian Cell (CHO)