

# Product Datasheet

## Human IL-27 ELISA Set 851.890.010

Product Name	Human IL-27 ELISA Set
Catalog Number	851.890.010
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	31.3-1000 pg/ml
Sensitivity	12 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-27
Sample Types	Supernatants, buffered solutions, serum, plasmas and other body fluids.