

## Product Datasheet

### Synuclein, alpha ( $\alpha$ -Syn) Antigen, Recombinant >85% LA662

<b>Product Name</b>	Synuclein, alpha ( $\alpha$ -Syn) Antigen, Recombinant >85%
<b>Catalog Number</b>	LA662
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Recombinant Antigen
<b>Storage</b>	$\leq -20^{\circ}\text{C}$ . Avoid repeated freezing and thawing.
<b>Purity</b>	>85% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 with 20% glycerol
<b>Matched Pairs</b>	HM1089, HM1092, HM1093 or HM1094
<b>Expression System</b>	E. coli
<b>Preservatives</b>	0.05% Sodium Azide