

Product Datasheet

Neurofilament L (NfL) Antigen, Recombinant >85% LA665

Product Name	Neurofilament L (NfL) Antigen, Recombinant >85%
Catalog Number	LA665
Tested Applications	LF, ELISA
Alternative Names	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS)
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	≤ -20°C. Avoid repeated freezing and thawing.
Purity	>85%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 20% glycerol
Molecular Weight	63 kDa
Matched Pairs	HM1095, HM1096, HM1097, HM1099, HM1100, HM1101
Expression System	Insect Cell
Preservatives	0.05% Sodium Azide