

Product Datasheet

Neurofilament L (NfL; 93-396) Antigen, Recombinant >85% LA667

Product Name	Neurofilament L (NfL; 93-396) Antigen, Recombinant >85%
Catalog Number	LA667
Tested Applications	ELISA, LF
Alternative Names	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS), Amyotrophic lateral sclerosis (ALS)
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	≤ -20°C. Avoid repeated freezing and thawing.
Purity	>85%
Buffer	Phosphate Buffered Saline, pH 7.2 with 20% glycerol
Molecular Weight	36.5 kDa
Matched Pairs	HM1095, HM1096, HM1097, HM1099, HM1100, HM1101, HM1184
Expression System	E. coli
Preservatives	0.05% Sodium Azide