

Product Datasheet

Neurofilament L (NfL; 93-396) Antigen, Recombinant >85% LA667

Product Name	Neurofilament L (NfL; 93-396) Antigen, Recombinant >85%
Catalog Number	LA667
Tested Applications	LF, ELISA
Alternative Names	NfL
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	≤ -20°C. Avoid repeated freezing and thawing.
Purity	>85%
Buffer	Phosphate Buffered Saline, pH 7.2 with 20% glycerol
Molecular Weight	36.5 kDa
Matched Pairs	HM1095, HM1096, HM1097, HM1099, HM1100, HM1101, HM1184
Expression System	E. coli
Preservatives	0.05% Sodium Azide