

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >85% LA780

Product Name Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen,
Recombinant >85%

Catalog Number LA780

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage 2-8°C

Epitope gp41

Purity >85%

Buffer 20 mM Carbonate-bicarbonate Buffer, pH 9.5 ± 0.2

Specificity gp41

Molecular Weight MW= 32 kDa

Matched Pairs LA782 (detection)

Expression System E. coli

Preservatives 0.1% Sodium Azide