

Product Datasheet

Human Immunodeficiency Virus 1 & 2 (HIV 1 & 2-gp36, gp41, gp120) Antigen, Recombinant >90% LA782

Product Name	Human Immunodeficiency Virus 1 & 2 (HIV 1 & 2-gp36, gp41, gp120) Antigen, Recombinant >90%
Catalog Number	LA782
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	2-8°C
Epitope	gp36, gp41, gp120
Purity	>90%
Buffer	20 mM Tris, pH 8.5 ± 0.2 with 5 mM EDTA
Specificity	gp36, gp41, gp120
Molecular Weight	MW= 43 kDa
Matched Pairs	LA780 or LA781 (capture)
Expression System	E. coli