

Product Datasheet

Hepatitis B Surface Antigen (HBsAg-subtype Ad), Native >95% LA787

Product Name Hepatitis B Surface Antigen (HBsAg-subtype Ad), Native >95%

Catalog Number LA787

Tested Applications ELISA

Brand EastCoast Bio

Form/Appearance Purified & 0.2 µm filtered Native Antigen

Storage 2-8°C

Note Tested negative for HIV 1/2 and HCV, however as with all biological material it should be handled as potentially infectious.

Purity >95%

Source Human Plasma

Buffer 25 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl, 5% sucrose

Preservatives 0.1% Sodium Azide