

Product Datasheet

HighScript Reverse Transcriptase #6601

Product Name	HighScript Reverse Transcriptase
Catalog Number	#6601
Description	HighScript Reverse Transcriptase is a modified version of MMLV reverse transcriptase with noticeable thermostability and high enzymatic activity. This enzyme is offered as a blend with an RNase inhibitor to prevent RNA degradation. HighScript Reverse Transcriptase together with enhanced buffer chemistry enables fast, efficient and accurate synthesis of the cDNA molecule.
Tested Applications	End-Point, Real-Time
Brand	myPOLS Biotec
Storage	-20°C
Content	10 000 U: 200 U/μl, 2 x 0.025 ml HighScript Reverse Transcriptase; 1 x 0.2 ml 5x HighScript buffer 40 000 U: 200 U/μl, 2 x 0.1 ml HighScript Reverse Transcriptase; 4 x 0.2 ml 5x HighScript buffer