

## Product Datasheet

### HighScript Reverse Transcriptase #6601

<b>Product Name</b>	HighScript Reverse Transcriptase
<b>Catalog Number</b>	#6601
<b>Description</b>	HighScript Reverse Transcriptase is a modified version of MMLV reverse transcriptase with noticeable thermostability and high enzymatic activity. This enzyme is offered as a blend with an RNase inhibitor to prevent RNA degradation. HighScript Reverse Transcriptase together with enhanced buffer chemistry enables fast, efficient and accurate synthesis of the cDNA molecule.
<b>Tested Applications</b>	End-Point, Real-Time
<b>Brand</b>	myPOLS Biotec
<b>Storage</b>	-20°C
<b>Content</b>	10 000 U: 200 U/μl, 2 x 0.025 ml HighScript Reverse Transcriptase; 1 x 0.2 ml 5x HighScript buffer  40 000 U: 200 U/μl, 2 x 0.1 ml HighScript Reverse Transcriptase; 4 x 0.2 ml 5x HighScript buffer