

Product Datasheet

MedixMDx RNase Inhibitor MX1105

Product Name	MedixMDx RNase Inhibitor
Catalog Number	MX1105
Description	MedixMDx RNase Inhibitor is a protein that blocks the activity of a variety of ribonucleases due to very fast non-covalent binding of the inhibitor to the active sites of enzymes. Medix RNase Inhibitor is resistant to oxidative stress due to its lack of cysteine residues and it remains active at 65°C for at least 30 minutes. This RNase inhibitor blocks a wide range of RNases, but does not inhibit RNases T1, T2, U1, U2, CL3, RNase I and H.
Brand	Medix Biochemica
Storage	-20°C
Content	2 500 U: 40 U/μl, 1 x 62.5 μl MedixMDx RNase Inhibitor