

## Product Datasheet

### qGreen Mix Separate ROX #2301

<b>Product Name</b>	qGreen Mix Separate ROX
<b>Catalog Number</b>	#2301
<b>Description</b>	qGreen Mix is a universal intercalating dye mix for robust, sensitive, and fast qPCR. qGreen Mix uses state-of-the-art technologies with an antibody-regulated hot-start Taq polymerase and intercalating dye for real-time PCR amplification. The optimized buffer chemistry and PCR enhancers and stabilizers enable rapid and sensitive qPCR.
<b>Tested Applications</b>	End-Point, Real-Time
<b>Brand</b>	Medix Biochemica
<b>Storage</b>	-20°C
<b>Content</b>	100 rxn: 1 x 1 ml qGreen Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive 500 rxn: 5 x 1 ml qGreen Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive