

Product Datasheet

Anti-Human CD25 Monoclonal Antibody, Biotin Conjugated Clone B-B10 852.004.010

Product Name	Anti-Human CD25 Monoclonal Antibody, Biotin Conjugated Clone B-B10
Catalog Number	852.004.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	PHA activated T cells
Brand	Diaclone
Form/Appearance	Lyophilised in phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. For longer storage freeze aliquots.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	PHA activated T cells
Specificity	Recognises the Interleukin 2 Receptor a chain (IL-2 Ra), a 55 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells