

# Product Datasheet

## MedixMDx qPCR lyo ready mix MX6101

Product Name	MedixMDx qPCR lyo ready mix
Catalog Number	MX6101
Description	MedixMDx qPCR Lyo ready mix is a universal probe mix that allows robust, sensitive, and fast qPCR for the detection of single or multiplex DNA targets. MedixMDx qPCR Lyo ready mix contains optimized excipients, buffer, PCR enhancers and an antibody-regulated hot-start Taq DNA polymerase and is ready to be lyophilized to produce stable reagents at room temperature. Upon addition of target specific primers/probes to the mastermix, the mixture can be lyophilized directly, without the need to add additional excipients. MedixMDx qPCR Lyo ready mix is compatible with several probes such as TaqMan® and Scorpions®.
Tested Applications	End-Point, Real-Time, Lyophilized
Brand	Medix Biochemica
Storage	-20°C
Content	500 rxn: 4 x 1.25 ml MedixMDx qPCR Lyo ready mix (2x) 10 000 rxn: 100 ml MedixMDx qPCR Lyo ready mix (2x)