

Product Datasheet

MedixMDx Anti-Taq Polymerase mAb (Clone MC1101) MX71031

Product Name	MedixMDx Anti-Taq Polymerase mAb (Clone MC1101)
Catalog Number	MX71031
Description	MedixMDx Anti-Taq Polymerase mAb is a highly purified monoclonal antibody that binds to DNA polymerase from Thermus aquaticus (Taq polymerase). Anti-Taq Polymerase mAb binds to and inactivates Taq DNA polymerase at room temperature. MedixMDx Anti-Taq Polymerase mAb dissociates from Taq DNA polymerase at elevated temperatures (higher than 50°C).
Brand	Medix Biochemica
Storage	+2°C - +8°C
Content	1 mg: 1 x 0.1 ml, 10 mg/ml MedixMDx Anti Taq Polymerase mAb 5 mg: 1 x 0.5 ml, 10 mg/ml MedixMDx Anti Taq Polymerase mAb