

Product Datasheet

Cortisol Antibody P01-92-94M-P

Product Name	Cortisol Antibody
Catalog Number	P01-92-94M-P
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Cortisol-3 CMO-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Corticosterone 45%, Progesterone 9%, Deoxycortisol <2, Dexamethazone <2%, Less than 0.01% to: Estrone (E1), Estriol (E3), Testosterone
Matched Pairs	P91-92-91H (detection)
Preservatives	0.05% Sodium Azide