

## Product Datasheet

### Methamphetamine Antibody P01-99-31R-IF

<b>Product Name</b>	Methamphetamine Antibody
<b>Catalog Number</b>	P01-99-31R-IF
<b>Description</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Methamphetamine-BgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	Methamphetamine-BgG linked at the Para Position off the Phenyl Group.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	Low cross reactivity to Amphetamine, Phenylpropanolamine
<b>Specificity</b>	Methamphetamine
<b>Matched Pairs</b>	P81-99-31A (capture) or P91-99-31H (detection)
<b>Preservatives</b>	0.05% Sodium Azide