

Product Datasheet

Secobarbital Antibody P01-99-41M-P

Product Name	Secobarbital Antibody
Catalog Number	P01-99-41M-P
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Secobarbital-KLH Linked off the 5 Position
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Buffer	0.05 M Phosphate, pH 7.15 and 0.85% Sodium Chloride
Specificity	Barbiturates
Matched Pairs	P81-99-41A (capture) or P91-99-41H (detection)
Preservatives	0.05% Sodium Azide