

Product Datasheet

Anti-h CRP 6402 SPTN-5 100145

Product Name	Anti-h CRP 6402 SPTN-5
Catalog Number	100145
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, FIA, IT
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Note	Note that CRP is a pentameric protein. The same antibody can thus be used both as a capture and a detection antibody.
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	$K_A = 6.2 \times 10^9 \text{ l/M}$; $K_D = 1.6 \times 10^{-10} \text{ M}$ (= 0.16 nM)
Associated Products	Recombinant CRP antigen, Lee Biosolutions 140-11R.
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.6-7.7
Cross Reactivity	Does not recognize serum amyloid-P complex
Specificity	Antibody recognizes human C-reactive protein (CRP), binding is Ca ²⁺ independent
Shelf Life	24 months
References	Käpyaho, K., Welin, M.-G., Tanner, P., Kärkkäinen, T. and Weber, T. (1989) Rapid determination of C-reactive protein by enzyme immunoassay using two monoclonal antibodies. Scand. J. Clin. Lab. Invest. 49:389-393