

Product Datasheet

Fentanyl Antibody P01-99-54R-IF

Product Name	Fentanyl Antibody
Catalog Number	P01-99-54R-IF
Description	Purified Polyclonal Antibody, IgG Fraction
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Fentanyl-BgG (linked off the aliphatic group, ethyl)
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Additional lots available, please inquire.
Clonality	Polyclonal
Buffer	0.05 M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Specificity	Fentanyl
Matched Pairs	P81-99-51A (capture) or P91-99-51H (detection)
Preservatives	0.05% Sodium Azide